



Adsorber Standing-Case



Product Code	ES	-	R	-	L
	Adsorber type case		Refillable cartridge		Size

Please note:

Temperature mustn't exceed 100°C during usage, or regeneration.

Dessicant	Filling weight [kg]
GIEBEL Xdry®	0.5
Silicagel	0.6
Active carbon	0.3
Molecular Sieve	0.4

Housing materials: SS, PVDF, PC, PA6

Total weight (kg): 0.437

NAME:

**Adsorber Standing-Case
ES-R_L**

Article no.:
18.0003.80

Status
freigabe

Maßstab

Format
A4

Blatt / von

1:1.5

1 / 1



GIEBEL Adsorber®
...setting standards in aeration drying!

rev.	Date	Name
	27.03.2024	Daniel Zinic

GIEBEL FilTec GmbH
www.giebel-adsorber.com
info@gf-dry.com