



Adsorber Standing-Case



| | | | | |
|--------------|--------------------|---|----------------------|----------|
| Product Code | ES | - | R | S |
| | Adsorber type case | | Refillable cartridge | Size |

Please note:

Temperature mustn't exceed 100°C during usage, or regeneration.

| Dessicant | Filling weight [kg] |
|-----------------|---------------------|
| GIEBEL Xdry® | 0.1 |
| Silicagel | 0.1 |
| Active carbon | 0.1 |
| Molecular Sieve | 0.1 |

Housing materials: SS, PVDF, PC, PA6

Total weight (kg): 0.138

NAME:

Adsorber Standing-Case
ES-R_S

Article no.:
18.0001.80

| | |
|-----------------|-------------------|
| Status freigabe | Format A4 |
| Maßstab 1:1 | Blatt / von 1 / 1 |



GIEBEL Adsorber®
...setting standards in aeration drying!

GIEBEL FilTec GmbH
www.giebel-adsorber.com
info@gf-dry.com